

Single Use Micro Cannula Needle For Dermal Filler ISO 13485

Our Product Introduction

for more products please visit us on hadermal-filler.com

Basic Information

- Place of Origin: China
- Brand Name: Blunt Micro Cannula
- Certification: CE
- Model Number: 18g-30g
- Minimum Order Quantity: 30pcs
- Price: 0.8usd/pcs
- Packaging Details: bag
- Delivery Time: 4-7work days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 13333box per week



Product Specification

- Product Name: Micro Cannula
- Material: Stainless Steel Blade
- Function: Micro Cannula For Dermal Filler Injection
- Application: Medical Disposable Blunt Tip Needle
- Usage: Single-use
- Package: 1pc/Pack,50pcs/box
- Type: Medical Disposable Blunt Tip Needle Micro Cannula
- Keyword: Blunt Tip Needle Micro Cannula
- Highlight: **Single use Micro Cannula Needle, microcannula for dermal filler ISO 13485, Single use microcannula for dermal filler**



More Images



Product Description

Medical Disposable Blunt Tip Needle Micro Cannula For Dermal Filler Injection

Micro cannula blunt tip needle for beauty injection

Specifications:

- 1) Superior quality and complies with international cosmetology and medical requirements.
- 2) Needle made of medical grade stainless steel.
- 3) sterile package
- 4) Disposable, for single use only.
- 5) Different size available, special design is available


What's The Procedure?


At the source of the filler injection using Microcannula, your Doctor creates a minute entry point in your skin using a precise, sharp needle. This is a vital step since the Microcannula has a blunt tip and is incapable of penetrating the skin by itself. The insertion location is minimal and does not cause substantial pain or bleeding. Then the Doctor skillfully inserts the Microcannula into the small opening and glides it under your skin to reach the sections of filler injection. The Microcannula then slides within the spaces connecting the tissue. It can be moved under the skin without any injuries to tissue. Patients typically do not feel much pain, if any, as the Doctor guides it under their skin.

Place of Origin	China
Disinfecting Type	EOS
Material	Stainless Steel
Safety standard	ISO 13485:2003
Application	Face, Neck, Chest and Body More



Item	Size	Needle Length	Product
Blunt Tip Needle (Micro Cannula)	18G	50mm	
		70mm	
		80mm	
	20G	100mm	
		38mm	
	22G	70mm	
		38mm	
	23G	50mm	
		38mm	
		40mm	
	25G	50mm	
		70mm	
		40mm	
	26G	50mm	
		70mm	
40mm			
27G	50mm		
	60mm		
	40mm		
30G	40mm		





shijiazhuang tiantian technology co.,ltd



008619931162089



info@tiantianbike.com



hadermal-filler.com

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA