



## Disposable 25mm 38mm Blunt Tip Microcannula Needle For Fillers

Our Product Introduction

### Basic Information

- Place of Origin: Korea
- Brand Name: Microcannula
- Certification: CE
- Model Number: Blunt Tip
- Minimum Order Quantity: 1 box
- Price: 40\$
- Packaging Details: box
- Delivery Time: 5-8days
- Payment Terms: T/T, Western Union, MoneyGram



### Product Specification

- Product Name: Microcannula
- Needle Size: 14g 16g 18g 19g 21g 22g 23g 25g 27g 29g 30g
- Length: 25mm 38mm 50mm 70mm 90mm 100mm
- Materials: Stainless Steel
- Function: Filler Injection
- Feature: Easy/Eco-friendly/Flexible
- Shelf Life: 5 Years
- Usage: Beauty Hospital Medical
- Highlight: Disposable blunt cannula needle, 25mm 38mm blunt cannula needle, blunt tip microcannula needle for fillers



### More Images



for more products please visit us on [hadermal-filler.com](http://hadermal-filler.com)

## Product Description

### Disposable Blunt Tip Cannula Blunt Microcannula Needle For Derma Filler Injection

#### What Is a Micro-Cannula?

A micro-cannula is a very thin, tiny, flexible hollow tube with a smooth round tip. The sharp end is gone! Instead of many needle sticks, we only need a very few to provide entry points for the cannula. That means fewer needle marks, less pain and far less probability of bruising.

#### Specification :

Product name	Microcannula Needle
Materials	Stainless steel
Size 18g	18gx50mm 18gx70mm 18gx100mm
Size 21g	21gx38mm 21gx50mm 21gx70mm
Size 23g	23gx38mm 23gx50mm 23gx70mm
Size 25g	25gx38mm 25gx50mm 25gx70mm
Size 27g	27gx38mm 27gx50mm 27gx70mm
Size 29g	29gx38mm 29gx50mm 29gx70mm
Size 30g	30gx25mm



18Gx50mm  
18Gx70mm  
21Gx50mm  
21Gx70mm  
22Gx50mm  
22Gx70mm  
23Gx50mm  
23Gx70mm  
25Gx38mm  
25Gx50mm  
27Gx38mm  
27Gx50mm  
30Gx25mm

#### What Are the Benefits of a Micro-Cannula over a Needle?

- Reduced pain
- Less bleeding
- Less bruising
- Less tissue trauma or swelling

Micro-cannula injection requires a single entry point for a large area instead of many as occurs with a traditional needle. Because the tip of the cannula is smooth and rounded, small vessels are simply pushed out of the way instead of being cut – so instead of pain you'll most frequently experience little more than pressure.

#### Why Use A Microcannula?

The decision to use a microcannula depends on the preference of the injector. Although needles are easy to use and can place product precisely where directed, they are more likely to cause bruising or transfer filler substance into a blood vessel. At depths where fillers are typically placed, most blood vessels cannot be seen. The utmost care will always be taken when injecting near areas that are more likely to bruise or with blood vessels we try to avoid. For instance, the tear trough region under the eye is more likely to bruise with a needle. Although not foolproof, microcannulas greatly reduce the risk of bruising and unintended injection into a vessel. Cannulas can also reduce trauma, discomfort and swelling.

#### So Why Use A Needle?


Since needles are more rigid and sharp, they can be inserted easily and placed just about anywhere with precision. They can be placed at all


depths of skin, under the skin and directly on top of bone. Cannulas can be injected at different levels of soft tissue, but not easily into the dermis of the skin. They are used when soft tissue implantation is desired beneath the dermis. A needle is sharp and easily penetrates the skin. It can be placed in the superficial, middle or deep layer of skin. Needles are also used for serial puncture methods of filler placement.

#### Do Fillers Come With Cannulas?

Fillers typically come packaged with at least two needles per syringe. These needles have sharp beveled ends that easily pierce the skin and soft tissue structures to precisely deliver the filler to its target. The needles can be substituted for others of different length and sizes. Cannulas are purchased separately, are more flexible and can be used instead of the supplied needles. Dr. Michael Lum uses a variety of cannula brands and sizes depending on the specific application.



 shijiazhuang tiantian technology co.,ltd

 008619931162089  [info@tiantianbike.com](mailto:info@tiantianbike.com)  [hadermal-filler.com](http://hadermal-filler.com)

NO 260 HUANAN ROAD SHIJIAZHUANG CITY HEBEI CHINA